


La nouvelle génération de bas anti-thrombose avec SEMELLE ANTIDÉRAPANTE (brevet européen)



Avantages

- Une meilleure élasticité et compression graduelle
- Utilisation sécurisante (prévention des chutes)
- Plus confortable pour le patient



6-8 mm/HG

12-14 mm/HG

18-20 mm/HG



Semelle
antidérapante



Talon
préformé



Ouverture au
devant du bas



Un bord en
silicone qualitatif



Extra
wide

Conseils de lavage

Lavage en machine Utiliser un filet de lavage. Peut être lavé à 90°C.
Eau de javel interdite. Séchage en machine interdit.
Ne pas repasser.

Lavage à la main Utiliser des détergents pour tissus délicats.
Rincer soigneusement, ne pas essorer.

La nouvelle génération des bas anti-thrombose avec SEMELLE ANTIDÉRAPANTE (brevet européen)





GUIDE DES TAILLES

ℓ = LONGUEUR
EXTRA WIDE

C = CIRCONFÉRENCE

	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE	XXX-LARGE
1 cB 18-22 cm	1 cB 22-26 cm	1 cB 26-29 cm	1 cB 29-32 cm	1 cB 32-35 cm	1 cB >35 cm	
2 cD 29-35 cm	2 cD 34-40 cm	2 cD 39-45 cm	2 cD 44-51 cm	2 cD 50-56 cm	2 cD >56 cm	
3 cG 43-54 cm	3 cG 48-60 cm	3 cG 54-66 cm	3 cG 60-72 cm	3 cG 67-80 cm	3 cG >80 cm	
3 cG 43-64 cm	3 cG 48-70 cm	3 cG 54-76 cm				

AD  1 cB	NORMAL 38-43 cm	H52 AS NS NORMAL SMALL	H52 AS NM NORMAL MEDIUM	H52 AS NL NORMAL LARGE	H52 AS NXL NORMAL X-LARGE	H52 AS NXXL NORMAL XX-LARGE
	LONG 43-48 cm 2 cD 1 cB	H52 AS LS LONG SMALL	H52 AS LM LONG MEDIUM	H52 AS LL LONG LARGE	H52 AS LXL LONG X-LARGE	-

AG  1 cB	SHORT 60-70 cm	H51 AS SS SHORT SMALL	H51 AS SM SHORT MEDIUM	H51 AS SL SHORT LARGE	H51 AS SXL SHORT X-LARGE	H51 AS SXXL SHORT XX-LARGE
	NORMAL 70-80 cm 3 cG 2 cD 1 cB	H51 AS NS H51 AS NSXW NORMAL SMALL	H51 AS NM H51 AS NMXW NORMAL MEDIUM	H51 AS NL H51 AS NLXW NORMAL LARGE	H51 AS NXL NORMAL X-LARGE	H51 AS NXXL NORMAL XX-LARGE
LONG 80-90 cm	H51 AS LS LONG SMALL	H51 AS LM LONG MEDIUM	H51 AS LL LONG LARGE	H51 AS LXL LONG X-LARGE	H51 AS LXXL LONG XX-LARGE	-